



CR
Carbon Rollator



For more joy of life:
Premium quality and fully equipped

The Saljol Carbon Rollator

Premium quality with many extras

The premium product impresses with its high quality and offers numerous features for your safety, such as over 30 reflectors, non-slip, skin-friendly handles, a practical stick holder, and top-quality brakes. With additional accessories, you can further optimise your rollator.

Safety you can feel – „hand-hugging“ Cork-TPE handles with armrest

The push handles are made from skin-friendly cork-TPE, which combines the sustainability, lightweight nature, and tactile feel of natural cork with the durability and resistance of plastics. They offer a secure grip thanks to their ergonomic wing shape, preventing your hands from slipping. The shock-absorbing and heat-retaining material reduces pain and promotes the mobility of your hands and fingers, especially for sensitive skin. The non-slip, quick-drying properties ensure secure comfort even in wet conditions. The Cork-TPE armrest makes resting your arms while seated very comfortable. As cork is a natural component of the handle, a typical dark patina may form due to the natural material's reaction with your palms. However, the handles are still easy to clean and remain hygienic.





Cork-coloured design elements guide you through the functions of Saljol Rollators

The design makes intuitive operation easier thanks to cork-coloured, brown components. All adjustable elements are immediately recognisable, allowing quick and secure adjustment to your needs. The well-thought-out colour scheme ensures clarity and maximum comfort in everyday use.



Visibility ensures safety

More than 30 reflectors and reflective strips on the rollator ensure that you remain highly visible in low-light conditions. Whether in traffic or during leisure activities, the reflective elements make you and your rollator more easily noticeable to other road users.



Optional one-hand braking system – brake with one or both hands for maximum control at no extra charge (on request)

The conversion to a one-hand simultaneous brake allows you to fully brake the rollator with just one hand. This braking system is particularly helpful for those experiencing age-related loss of strength or unilateral limitations, such as after a stroke, with Parkinson's, or fractures, and enables you to operate both brakes with minimal effort. The necessary longer cable for the brake system will be sent to you free of charge upon request.



Compact and foldable – ideal for easy transport

The rollator can be folded to a width of just 21 cm with a single movement. When folded, it can be securely locked with a hook, making it easy to load and store compactly and stably. Thanks to its lightweight design and narrow dimensions, it can be easily transported in a car or while travelling.



Highest quality in every detail – you'll feel it!

The smooth-running hand brakes allow for effortless braking and are easy to operate, even with smaller hands. The ergonomic shape and cork-TPE covering provide a secure and comfortable grip. When you engage the parking brake, a clear clicking sound lets you know that the rollator is securely in place – ideal for sitting down or filling the bag.



Built-in safety: grinding brake for better control when going downhill or when more support is needed

Safety and control while walking are your top priorities? The grinding brake on the CR provides exactly that. The rolling resistance of the rear wheels can be individually adjusted, gently braking without locking (ABS), and preventing your rollator from rolling faster than you want – especially helpful on slopes. The brake supports you optimally in cases of walking instability, can counteract incorrect directional drift, and thus reduces the risk of falling. Once properly adjusted, you can activate or deactivate it with your foot.



Very light and stable with a total weight of only 5.6 kg

The Carbon Rollator is exceptionally light and manoeuvrable, weighing just 5.6 kg and with a width of 61 cm. The weight has been carefully chosen as a good compromise between driving stability and ease of handling. Too light can quickly become unstable, while too heavy can be an obstacle when loading or carrying.



Secure hold for your walking stick – always within reach

The universal cane holder adapts flexibly to different mobility aids and walking sticks. With a unique locking system, the stick is securely fixed and always within reach. A clicking sound lets you know that the stick is safely attached.



Your support wherever you go – the lightweight and foldable walking stick

The walking stick is incredibly lightweight at just 278 g, yet remarkably durable. Made from high-quality carbon fibre, it impresses with its top quality and elegant look. It is also height-adjustable, foldable and always stable – the perfect combination of functionality and elegance. It is available directly at www.saljol.de/en, through our service team, or at a medical supply store near you.



The right product for you – available for every need

Saljol Premium Rollators are available in two sizes, elegant colours, and additionally with a wider seat (52 cm) and higher weight capacity (180 kg). This gives you the option to choose the solution that best suits your needs and personal taste. Whether for more comfort, functionality, or style – you will find the perfect combination to meet your requirements.



Image shows AR62W





German Engineering from Bavaria

German Design



Great luxury that remains affordable

Saljol Rollators offer the perfect mix of functionality, quality, and a touch of luxury. German engineering combines high-quality features with comfort to make your everyday life easier – all at a fair overall price.



Comprehensive support – for your mobility and safety

Our physiotherapist, Katharina, supports you in all areas where expert physiotherapeutic knowledge is required. Our helpful customer service is always at your side – whether in the medical supply store or directly with us. If you have any questions or uncertainties regarding your Saljol Carbon Rollator, we provide the support you need to get the most out of your rollator.



Training with your rollator

To avoid medical side effects when using our rollators, we invented Rollastic. It's a tailored exercise programme for rollator users that, when performed regularly, can help maintain mobility and strength. Special videos are provided with the rollator. For your safety in road traffic, we developed the safety training „Rollator Mobility“. Helpful tips and tricks can be found at www.rollatormobility.com.



Marietta: „I enjoy showing off my Rollator!“

„I suffer from HMSN/CMT, a neuromuscular disease. That's why I like the ergonomically shaped handles of the Carbon Rollator, as they give me more stability. The wide support surface allows me to safely lean on the rollator, especially when I experience pain in my legs on hard surfaces. Unlike plastic handles, rain doesn't affect these handles. You don't slip, and you don't feel like your hands are wet. The safety is still there, and they don't get cold quickly. Also, my carbon rollator is my Porsche – I can go out in style with it. I have no reservations about that.“



Carbon Rollator	CR54	CR62	CR62W
Recommended body height	145–170 cm	160–200 cm	160–200 cm
Maximum user weight	150 kg	150 kg	180 kg
Total weight without accessories	5.6 kg	5.6 kg	5.6 kg
Maximum load of bag	5 kg	5 kg	5 kg
Width open/closed	61/21 cm	61/21 cm	68/21 cm
Seat width	45 cm	45 cm	52 cm
Seat height	54 cm	62 cm	62 cm
Handle height (adjustable)	76–88 cm	81–103 cm	81–100 cm

Carbon Rollator accessories



Back belt



LED light



Bell



Stick holder



Walking stick



Firm seat



Travel bag



Wooden tray
with felt pad



Umbrella
caramel



Umbrella
black



JOST-bag
black



Brake cable for
the conversion to
one-hand brake

Frame colour



British Racing
Green



Midnight Blue



Copper Red



Pure Gold