

Saljol Carbon Rollator

The safe and
elegant companion



Highlights

Skin-friendly handles

cork-TPE provides a natural feeling, is shock-absorbing and gentle on arms and wrists

Bag can be removed with one hand

more safety – hold on to the rollator handle with one hand, unhook the bag with the other hand

Shock absorbing all-round tyres

ideal for town and country – special tyres for more grip, comfort and safety on any ground



Armrest made of cork-TPE

sit comfortably without pressure points on the forearm


Grinding brakes

prevent the rollator from rolling away when walking downhill, upright walk in the rollator is supported and directional twist is compensated when walking

Reflectors on all sides

good visibility thanks to over 30 reflectors on frame, tyres, handles, tipping aid and bag

Easy
transport



Sit stably
and
comfortably



Go downhill
safely with
the grinding
brake



Technical data

Carbon Rollator	CR54	CR62
Recommended for body height	<170 cm body height	>160 cm body height
Maximum user weight	150 kg	150 kg
Weight without bag	5.6 kg	5.6 kg
Weight with bag	6.0 kg	6.0 kg
Maximum load of the bag	5 kg	5 kg
Width open	61 cm	61 cm
Width closed	21 cm	21 cm
Seat width	45 cm	45 cm
Seat height	54 cm	62 cm
Maximum handle height	88 cm	103 cm
Minimum handle height	76 cm	81 cm

Frame colour Carbon Rollator



British Racing Green



Midnight Blue



Star Silver



Matt Black



Conversion set for one-hand simultaneous brake

Accessories



Transport bag



Back belt



Designer bag from JOST

leopard print or black




Bag



Firm seat



Real wood tray



Parasol/umbrella



Additional cane holder



Bell



LED lamp