



# CR

Carbon Rollator



For more joy of life – premium quality and fully equipped

**Includes back belt, bag, stick holder, LED light and bell**

Safety you can feel –  
'hand-friendly' cork TPE  
handles with armrest

Effortless braking with  
one hand – convertible  
to a one-hand brake  
system

Grinding brake for greater  
control when riding  
downhill and if you feel  
unsteady when walking



Weights only  
5.6 kg

Compact and  
foldable – easy  
to transport

Visible even in  
the dark thanks  
to over  
30 reflectors

Highest quality down to the  
smallest detail – you can feel it!

Carbon Rollator	CR54	CR62	CR62W*
Recommended body height	145–170 cm	160–200 cm	160–200 cm
Maximum user weight	150 kg	150 kg	180 kg
Total weight without accessories	5.6 kg	5.6 kg	5.6 kg
Maximum load of bag	5 kg	5 kg	5 kg
Width open/closed	61/21 cm	61/21 cm	68/21 cm
Seat width	45 cm	45 cm	52 cm
Seat height	54 cm	62 cm	62 cm
Handle height (adjustable)	76–88 cm	81–103 cm	81–100 cm

\*Only available in British Racing Green

### Carbon Rollator accessories



Back belt



LED light



Bell



Stick holder



Firm seat



Travel bag



Wooden tray  
with felt pad



JOST-bag  
black



Umbrella  
caramel



Umbrella  
black



Brake cable for  
the conversion to  
one-hand brake

### Frame colour



British Racing  
Green



Midnight Blue



Copper Red



Pure Gold

Products may differ from illustrations | 07/2025