



CR

Carbon Rollator



For more joy of life – premium quality and fully equipped

Includes back belt, bag, stick holder, LED light and bell

Safety you can feel –
'hand-friendly' cork TPE
handles with armrest

Effortless braking with
one hand – convertible
to a one-hand brake
system

Sliding brake for greater
control when riding
downhill and when
unsure of the gear



Weights only
5.6 kg

Compact and
foldable – easy
to transport

Visible even in
the dark thanks
to over
30 reflectors

Highest quality down to the
smallest detail – you can feel it!

Carbon Rollator	CR54	CR62	CR62W*
Recommended body height	145–170 cm	160–200 cm	160–200 cm
Maximum user weight	150 kg	150 kg	180 kg
Total weight without accessories	5.6 kg	5.6 kg	5.6 kg
Maximum load of bag	5 kg	5 kg	5 kg
Width open/closed	61/21 cm	61/21 cm	68/21 cm
Seat width	45 cm	45 cm	52 cm
Seat height	54 cm	62 cm	62 cm
Handle height (adjustable)	76–88 cm	81–103 cm	81–100 cm

*Only available in British Racing Green

Carbon Rollator accessories



Back belt



LED light



Bell



Stick holder



Walking stick



Firm seat



Travel bag



Wooden tray
with felt pad



Umbrella
caramel



Umbrella
black



JOST-bag
black



Brake cable for
the conversion to
one-hand brake

Frame colour



British Racing
Green



Midnight Blue



Copper Red



Pure Gold

Products may differ from illustrations | 07/2025