

Rise & Recline Chairs

Sit and lie down comfortably –
take a seat and stand up safely



CE

Rise & Recline Chair

Depending on the version, you can enjoy these comforts:

- Lumbar support for additional stability when sitting
- magnetic remote control easy to attach directly to the chair
- get into a sitting, reclining or standing position with the remote control
- stand-up support in the armrest
- adjustable headrest
- adjustable seat depth
- soft or firm upholstery – You have the choice



Fall prevention

- ✓ magnetic remote control prevents dangerous falls
- ✓ a rechargeable battery as an accessory avoids tripping hazards due to cables lying around

Club1: stand up easily, easy to use



Club2: sleep and relax in one armchair



Club2 Soft: Armchair with special comfort



Club2 Mobil: the functional armchair




Club3: relief for the legs due to heart balance position



Club1: simple handling with only two buttons

Maximum
lying position



Magnetic remote control

adheres to the side and on top of the armrest and to the side of the headrest – no more bending down for the remote control

Lumbar support

additional stability when sitting



Simple handling

get into all positions with only two buttons

Seat depth adjustment

individually adjustable to your leg length

Spring core upholstery

sit comfortably


Colour variants



Microfibre
beige



Woven
fabric grey



- only two buttons
- simplest handling

Magnetic remote control



Stand up safely

Relax comfortably



*GGT stands for Gesellschaft für Gerontotechnik (Society for Gerontological Technology); the test seals of the neutral, German, non-profit competence carrier confirm that products are particularly recommendable in terms of function, user-friendliness, maintenance and care.

Club2: sleep and relax in one armchair

**Maximum
lying position**



Adjustable headrest

more comfort due to
tiltable headrest

Magnetic remote control

adheres to the side and on
top of the armrest and to
the side of the headrest –
no more bending down
for the remote control

Lumbar support

additional stability
when sitting

Easy change of position

easy to move in the room
thanks to castors

Stepless adjustment of the backrest

armchair stops automatically
as soon as the remote
button is released

Seat depth adjustment

individually adjustable to
your leg length

Day bed function

take a nap without having
to change to the bed, side
panel prevents falling out

Separate control of the footrest

put your feet up,
while the back
remains straight


Colour variants



Woven fabric grey &
Imitation leather dark brown



Microfibre beige &
Imitation leather dark brown



- turns into a day bed at the
touch of a button
- backrest and footrest
separately adjustable

Stand up and sit down
safely



Separate control of
backrest and footrest

Lying horizontally
like in a bed



Club2 Soft: sleeping chair with extra comfort

Maximum
lying position



Adjustable headrest

more comfort due to
tiltable headrest

Magnetic remote control

adheres to the side and on
top of the armrest and to
the side of the headrest –
no more bending down
for the remote control

Lumbar support

additional stability
when sitting

Extra soft due to cushions

removable and adjustable
cushions – ideal for people
with sensitive skin

Infinitely variable adjustment of the backrest

armchair stops automatically
as soon as the remote
button is released

Day bed function

take a nap without having
to change to the bed, side
panel prevents falling out


Separate control of the footrest

put your feet up,
while the back
remains straight

Seat depth adjustment

individually adjustable to
your leg length

Colour variants




Microfibre
beige



Woven
fabric grey

Easy change of position

easy to move in the room
thanks to castors



- turns into a day bed at the touch of a button
- backrest and footrest separately adjustable

Relax in the Club2 Soft



Soft and comfortable
thanks to extra padding


Four castors for easy moving



Club2 Mobil: the functional armchair

All functions of Club2 included

Maximum
lying position



Push handle

facilitates the transport
of persons

Battery as accessory

more mobility in the
living area, no tripping
hazard due to power
cable lying around

Extra large, stable castors

facilitate the transport
of persons


Hygienic

can be disinfected by wiping and
is water-repellent

Removable armrests

for the „sideways“ transfer
of persons from the bed/
wheelchair into the chair

Colour variants



Imitation leather dark brown




Woven fabric grey &
Imitation leather dark brown

Wipe disinfection possible!

The Club2 Mobil in imitation
leather version can be
disinfected by wiping

Fire resistant according
to BS standard 5852 to
foam and cover fabric



- programmable favourite position
- illuminated buttons
- backrest and footrest separately adjustable

With push handle for
easy transport



Easy transfer from
wheelchair into the armchair



Removable armrests



Club3: Riser and Recliner Armchair for maximum relaxation

**Maximum
lying position**



Adjustable headrest

more comfort due to
tiltable headrest

Lumbar support

additional stability
when sitting

Heart balance position

the position with the legs
above the heart relieves
the pressure on your legs

4 speed-dial buttons

for quick setting of
the desired position

Little space required


installation area is
only 150 x 75 cm –
ideal for small flats and
little space

Magnetic remote control

adheres to the side and on
top of the armrest and to
the side of the headrest –
no more bending down
for the remote control

Stepless seat tilt adjustment

tiltable up to 108°



- programmable favourite position
- illuminated buttons
- backrest and footrest separately adjustable

Colour variants




Woven fabric grey &
Imitation leather dark brown



Microfibre beige &
Imitation leather dark brown

Relaxation due to
pressure relief



Stand-up function

Relaxation in the
heart balance position



Accessories



Protective seat covers

- contains 1 cover for footrest and 1 cover each for armrest left and right
- washable up to 40°C
- available in the colours
Microfibre beige
Woven fabric grey
Imitation leather dark brown

not available for Club2 Soft



Incontinence pad set

- contains 2 cushions for the seat
- no slipping and easy to replace thanks to velcro system
- washable up to 40°C
- protects with an integrated 100% moisture barrier made of PU membrane and a cotton layer that absorbs up to 500 ml of liquid
- available in the colours
Microfibre beige
Woven fabric grey



Lumbar cushion

- supports the lumbar spine
- can also be used for soft positioning of the feet and heels
- available in the colours
Microfibre beige
Woven fabric grey
Imitation leather dark brown

Available accessories may vary depending on the country.

Accessories



Head cushion

- supports the cervical spine
- the head cushion can be adapted individually to the neck
- available in the colours
Microfibre beige
Woven fabric grey
Imitation leather dark brown



Battery

- for wireless use of the chair
- avoids tripping hazards due to cables lying around
- more mobility with the living area

Head cushion accessory



Available accessories may vary depending on the country.

The swivel table for all Club-chairs


- the perfect support: holds up to 20 kg on the table surface
- very stable and supportive: you can lean on the tube with up to 100 kg
- when using the stand-up function, the table automatically turns to the side and thus cannot damage the chair
- with the help of two rings at the foot of the table, the height of the table can be adjusted
- is firmly screwed to the armchair, can be mounted on the left or right side
- 270° swivelling, swivelling radius 70 cm
- can be disinfected by wiping
- available in bright and dark wood look



Swivel table bright



Swivel table dark
























turns away automatically
when using the stand-up
function



Large table top

Firmly screwed to
the armchair



Find the armchair that suits you:		Colour variants		Inclination backrest	Seat width (cm)	Seat height (cm)	Seat depth (cm)	Height (cm)	Chair height (cm)	Chair width (cm)	Chair length in rest position (cm)	Maximum load (kg)	Floor space (cm)
Rest and stand-up chair	All Club rest and stand-up chairs can be retrofitted with a rechargeable battery and have the following functions: <ul style="list-style-type: none">• magnetic remote control• haptic support point in the armrest on decorative stitching• lumbar support• stand-up function												
Club1	Simplest operation with only two buttons. Additional functions: <ul style="list-style-type: none">• super comfortable thanks to spring core• seat depth adjustment to adapt the seat to your body measurements	 Microfibre beige	 Woven fabric grey	133°	55	48	56/61	160-180	112	75	170	150	180x75
Club2	Ideal for sleeping and relaxing. Additional functions: <ul style="list-style-type: none">• day bed function: horizontal reclining possible• adjustable headrest• rollable on 4 wheels• seat depth adjustment to adapt the seat to your body dimensions• separate control of the footrest	 Woven fabric grey & Imitation leather dark brown	 Microfibre beige & Imitation leather dark brown	166°	55	46	47/51	150-190	121	75	194	150	200x75
Club2 Soft	Even more comfortable with extra soft padding, ideal for thin and sensitive skin. Additional functions: <ul style="list-style-type: none">• day bed function: horizontal lying possible• adjustable headrest• removable upholstery• rollable on 4 wheels• seat depth adjustment to adapt the seat to your body dimensions• separate control of the footrest	 Microfibre beige	 Woven fabric grey	166°	49	50	47/51	150-190	123	75	194	150	200x75
Club2 Mobil	Large wheels and push handle for transport. Additional functions: <ul style="list-style-type: none">• day bed function: horizontal reclining possible• can be rolled on 4 wheels and is approved for the transport of persons• can be wipe-disinfected and is water-repellent• removable armrests for lateral transfer of persons from bed/wheelchair to chair• seat depth adjustment to adapt the seat to your body measurements <ul style="list-style-type: none">• Adjustable headrest• separate control of the footrest	 Imitation leather dark brown	 Woven fabric grey & Imitation leather dark brown	166°	55	50	47/51	150-190	128	75	190	150	200x75
Club3	Armchair for maximum relaxation. Additional functions: <ul style="list-style-type: none">• relief for the legs due to heart balance position• adjustable headrest• programmable favourite position• illuminated magnetic remote control• separate control of the footrest	 Woven fabric grey & Imitation leather dark brown	 Microfibre beige & Imitation leather dark brown	108°	53	50	54	150-190	119	75	151	150	150x75

