

Saljol Allround Rollator

The all-rounder
for city and countryside



Highlights

Convertible brake

standard two-hand brake
or one-hand simultaneous
brake

Aluminium frame with new design

lightweight frame made of
aluminium in Pearl Grey

All-terrain tyres

extra-large 250 mm
front wheels for
more driving comfort

Light and robust

total weight only 6.4 kg



Armrest made of cork-TPE

sit comfortably without
pressure points on the
forearm

Lightweight net bag

loadable with 5 kg

Reflectors

good visibility on all sides
due to over 30 reflectors
on frame, tyres, handles,
tipping aid and bag

Improved tipping aid

anti-slip pads on the
tipping aid facilitate
tilting of the rollator

Cork-TPE armrest



Extra-large front wheels

Convertible to
one-hand simultaneous brake



Cork-Highlights

Skin-friendly cork TPE on the push handles

- shock absorbing and heat absorbing
- UV-resistant and abrasion-proof
- non-slip even with sweaty hands



Functional elements are recognizable by cork



Technical data

Allround Rollator	
Recommended body height	> 160 cm
Maximum user weight	150 kg
Total weight without net bag	6.4 kg
Maximum net bag load	5 kg
Width (open)	61 cm
Width (closed)	21 cm
Seat width (between push handles)	46 cm
Seat height	62 cm
handle height (adjustable)	81–103 cm

Accessories



Travel bag



Back belt



Jost bag leopard print



Jost bag black



Standard bag



Bag black/brown



Real wood tray




Firm seat



LED lamp



Bell



Parasol/umbrella



Additional cane holder

Products may differ from illustrations, version 08/2022