



AR

Allround Rollator



For more joy of life – premium quality and fully equipped

Includes back belt, net bag, stick holder, LED light and bell

Safety you can feel –
'hand-friendly' cork TPE
handles with armrest

Effortless braking with
one hand – convertible
to a one-hand brake
system

Tyres for all terrains
– shock-absorbing,
puncture-proof and
extra-large front wheels



Weights only
6.4 kg

Compact and
foldable – easy
to transport

Visible even in
the dark thanks
to over
30 reflectors

Highest quality down to the
smallest detail – you can feel it!

Allround Rollator	AR54	AR62	AR62W*
Recommended body height	145–170 cm	160–200 cm	160–200 cm
Maximum user weight	150 kg	150 kg	180 kg
Total weight without accessories	6.3 kg	6.4 kg	6.4 kg
Maximum load of net bag	5 kg	5 kg	5 kg
Width open/closed	61/21 cm	61/21 cm	68/21 cm
Seat width	46 cm	46 cm	52 cm
Seat height	54 cm	62 cm	62 cm
Handle height (adjustable)	76–88 cm	81–103 cm	81–100 cm

*Only available in Barolo Red

Allround Rollator accessories



Back belt



LED light



Bell



Stick holder



Grinding brake

Illustration shows CR



Travel bag



Firm seat



Wooden tray
with felt pad



Umbrella
caramel



Umbrella
black



Walking stick



Brake cable for
the conversion to
one-hand brake

Frame colour



Pearl Grey



Magic Black



Forest Green



Barolo Red

Products may differ from illustrations | 07/2025