

# Saljol Shower Stool

Turns safely  
on the spot



Product innovation - DHK59:  
easier to stand up thanks  
to wedge seat



## Highlights



### **More safety**

when showering after operations or  
Standing/stand-up insecurity

### **Heat-reflecting seat**

comfortable feeling when sitting

### **Maintaining privacy**

Self-determined and safe washing

### **New: shower stool also available with wedge seat**

supports the shower user when standing up

### **Safe 360° turn on the spot**

effortless turning of the stool makes it easier  
to reach the fittings

### **Open front frame**

offers more tilt resistance and facilitates entry

### **Anti-slip wheels**

arrangement of wheels made of non-slip TPE  
material prevents the stool from slipping away and  
provides safety when stepping into the shower tray

## Technical data

Saljol Shower Stool	DH51	DH55	DHK59
Recommended for body height	1.50-1.70 m	>1.70 m	>1.50 m
Maximum user weight	150 kg	150 kg	150 kg
Hygiene cut-out, setting 1/2/3	5/8/11 cm	5/8/11 cm	5/8/11 cm
Total width from handle to handle	52 cm	52 cm	52 cm
Seat height	51 cm	55 cm	59 cm (middle)
Wheelbase	32 cm	32 cm	32 cm
Space required for turning circle	48 cm	48 cm	48 cm
Total weight Velvet Grey/Diamond Silver	5.8/4.4 kg	6.0/4.6 kg	6.0/4.7 kg

## Frame colour



Velvet Grey



Diamond Silver

Version with wedge seat



Triple adjustable hygiene cut-out



Seat reflects body heat



Products may differ from illustrations, version 08/2022

11/2022