

# Saljol Shower Stool

Turns safely  
on the spot



Shower Stool with  
wedge seat – DHK59:  
easier to stand up

CE

 **SALJOL**  
Joy of life!

## Highlights



### **More safety**

when showering after operations or  
standing/stand-up insecurity

### **Heat-reflecting seat**

comfortable feeling when sitting

### **Maintaining privacy**

self-determined and safe washing

### **Shower stool also available with wedge seat**

supports the shower user when standing up

### **Safe 360° turn on the spot**

effortless turning of the stool makes it easier  
to reach the fittings

### **Open front frame**

offers more tilt resistance and facilitates entry

## Technical data

| Saljol Shower Stool                     | DH51        | DH55       | DHK59          |
|---|-------------|------------|----------------|
| Recommended for body height             | 1.50-1.70 m | >1.70 m    | >1.50 m        |
| Maximum user weight                     | 150 kg      | 150 kg     | 150 kg         |
| Hygiene cut-out, setting 1/2/3          | 5/8/11 cm   | 5/8/11 cm  | 5/8/11 cm      |
| Total width from handle to handle       | 52 cm       | 52 cm      | 52 cm          |
| Seat height                             | 51 cm       | 55 cm      | 59 cm (middle) |
| Wheelbase                               | 32 cm       | 32 cm      | 32 cm          |
| Space required for turning circle       | 48 cm       | 48 cm      | 48 cm          |
| Total weight Velvet Grey/Diamond Silver | 5.8/4.4 kg  | 6.0/4.6 kg | 6.0/4.7 kg     |

## Frame colour



Version with wedge seat



Triple adjustable hygiene cut-out



Seat reflects body heat



Products may differ from illustrations

01/2025