

Saljol Comfort Mattress

Sleep relaxed

with
KUBIVENT
inside

CE



Highlights

Choice of removable cover

Terrycloth or incontinence cover, washable at 60° C

Three foams for high comfort

Optimal support of the spine, good weight distribution and promotion of air circulation

Suitable for care beds

mattress can be fitted in most care beds

Available in different sizes

in standard and special sizes obtainable



Cosy and soft



Three foam layers for extra high comfort

Skin friendly cover




Sizes

The Saljol Comfort Mattress is available in almost any size:

Standard dimensions:	90 x 200 cm (Single Size)	Special dimensions:	90 x 190 cm	90 x 220 cm
	140 x 200 cm (Double Size)		100 x 200 cm	120 x 200 cm
	160 x 200 cm (Queen Size)		120 x 220 cm	and more

Mattress thickness 15 cm. Additional sizes on request. Available within 10 working days.

Structure and indication



- 1 2 cm soft foam top layer**
(volume weight 35 kg/m³, compression hardness 1.8 kPa)
- 2 4 cm cold foam nubbed pad**
(volume weight 45 kg/m³, compression hardness 3.5 kPa)
- 3 9 cm cold foam bottom layer**
(volume weight 40 kg/m³, compression hardness 3.0 kPa)

Recommended for

- normal sleepers without particular sleep problems
- back pain
- neck pain
- headaches due to tension

Versions

Incontinence cover
does not allow moisture to pass through

Terrycloth cover
particularly soft and skin friendly



Made in Germany

Products may differ from illustrations

12/2022