

Saljol Carbon Walking Stick

Elegant, ergonomic and foldable

CE



Highlights

Stylish and lightweight helper

chic look thanks to glossy carbon and low weight of only 278 grams

Ergonomic and safe

Grip handles with cork TPE provide a natural feel, are shock-absorbing and gentle on arms and wrists, pole length adjustable from 83.5 to 93.5 cm

Robust support aid

Made of sturdy carbon and loadable up to 100 kg



Two grip areas

upper handle for use as a support pole, lower handle for support when walking

Integrated safety wrist strap

the pole hangs securely on the wrist

Easy to transport thanks to folding mechanism

Simply pull apart and fold, pack size when folded: 43 x 12 x 3 cm

Easy to fold



Elegant and practical like all Saljol products

Useful accessories for free hands - see back



Saljol magnetic cane holder

No need to bend down again

CE



For on the road: Held on a stick, Taler in your pocket



Do you use a walking stick or crutches? And sometimes the walking aid falls to the ground? Bending down is strenuous and wobbly. You could fall. This is where our magnetic cane holder Held and the magnetic Taler help: you can attach your walking aid to your trousers, for example.



Held
black, yellow, orange




Taler
black, yellow, orange




Held/Taler set
black, yellow, orange

Metal adapter for at home



The Held set consists of a magnetic stick holder and three self-adhesive metal adapters: For your favourite place at the table, as a wall mount and for universal use.

The Held complete set consists of a magnetic cane holder, a Taler for on the go and three self-adhesive metal adapters.



Held set for indoors



Held complete set

Products may differ from illustrations, version 11/2022

12/2022