

Saljol Carbon Walking Stick

Elegant, ergonomic and foldable

CE



Highlights

Stylish and lightweight helper

chic look thanks to glossy carbon and low weight of only 278 grams

Ergonomic and safe

Grip handles with cork TPE provide a natural feel, are shock-absorbing and gentle on arms and wrists, pole length adjustable from 83.5 to 93.5 cm

Robust support aid

Made of sturdy carbon and loadable up to 100 kg



Two grip areas

upper handle for use as a support pole, lower handle for support when walking

Integrated safety wrist strap

the pole hangs securely on the wrist

Easy to transport thanks to folding mechanism

Simply pull apart and fold, pack size when folded: 43 x 12 x 3 cm



Easy to fold



Easy to stow away



Elegant and practical like all Saljol products