

# Saljol Indoor Rollator

Stable helper  
for your living area

MD CE



## Highlights

### Move safely in the home

invisible brake cables enable driving without getting caught on furniture or door handles

### Stable handrail

lean on it safely without the rollator tipping over

### 360° rotatable wheels

turn on the spot, even in confined spaces, such as the guest toilet



### Innovative braking system

rollator offers two braking options: one- or two-handed at any position of the brake bar or two-handed with the brake handles

### Convenient wooden seat

move around the flat while sitting and save energy


### Narrow design

easy to manoeuvre in tight spaces, also fits through narrow door frames


Replaces toilet handhold



Fits through narrow spaces



Enables tripping while sitting



# Technical data

Indoor Rollator	WR52	WR55	WR58	WRXL
Recommended for body height	145-165 cm	150-180 cm	170-200 cm	170-210 cm
Maximum user weight	150 kg	150 kg	150 kg	200 kg
Handle height	80 cm	86 cm	92 cm	95 cm
Seat height	52 cm	55 cm	58 cm	58 cm
Width at the bottom (wheel)	56 cm	56 cm	56 cm	68 cm
Width at the top	51 cm	51 cm	51 cm	65 cm
Distance between handles	43 cm	43 cm	43 cm	60 cm
Weight Velvet Grey/Ivory White without accessories	10.4 kg	10.8 kg	11.2 kg	13.4 kg
Weight Diamond Silver without accessories	9.3 kg	9.5 kg	9.7 kg	
Compatible with accessories	✓	✓	✓	

## Frame colour



Ivory White

Velvet Grey

Diamond Silver

All models have a removable brake bar in the respective frame colour.

The Indoor Rollator in Diamond Silver is made of electro-polished stainless steel and is also supplied with a soft seat and back cushion.

## Accessories



Tray



Jost bags



Jost bags available in leatherette beige, transparent or felt grey



Back roll



Upholstered seat



Cabinet

WRXL indoor rollator available in Velvet Grey

Also available in Ivory White on request