



# RA

Riser Armchair





Safe and sound:  
your perfect spot for relaxation

## The Saljol Riser Armchair

### Premium quality with lots of features

This premium product stands out with its high quality and offers a range of safety-enhancing features, including a rise and recline function, a magnetic remote control, a side pocket, and a USB port. With additional accessories, you can further customise your chair to suit your needs.

### Smart and elegant design for your well-being

The Saljol Riser Armchair combines a sleek, modern look with maximum comfort. Its thoughtfully designed structure creates an inviting atmosphere while providing a sense of security and cosiness. The distinctive „wings“ offer additional privacy, allowing you to relax even more deeply. A true feel-good piece of furniture that perfectly combines style and function!







### Quality that speaks for itself

This chair perfectly combines comfort and premium materials. High-quality comfort foam provides ideal support and a pleasantly soft feel – perfect for both sitting and reclining in total relaxation.

The upholstery is easy to clean, water-repellent, fire-resistant, urine-resistant, and machine washable at 60 °C – ensuring maximum hygiene and safety. All materials comply with strict EU standards (REACH, RoHS).



### Adjustable seat depth – perfectly customised

The riser armchair adapts effortlessly to your leg and thigh length with its two-stage adjustable seat depth. Enjoy maximum comfort, a healthy sitting posture, and the perfect relaxation experience – whether you're unwinding or getting up with ease.

### Two colours – perfect for your home

Platinum Grey and Graphite Grey – the elegant chair is available in these two stylish shades. Like a chameleon, it blends seamlessly into any living space. Choose your upholstery set in either of these colours to suit your taste and turn your chair into an eye-catcher in your personal haven of comfort.



### Choose the upholstery that suits you best

Since the base of the riser armchairs is essentially the same, you have the freedom to select your preferred upholstery: Classic upholstery – elegant and supportive for a firmer, more stable sitting experience or soft upholstery – plush and cosy for maximum comfort.

Both models offer premium comfort and an exceptional seating and reclining experience as the seat area and lumbar section of the backrest have been specially reinforced. Whether sitting or reclining, your back experiences noticeable relief – especially if you spend a long time in the armchair. By varying the cushion types, you can easily adjust the seating position to suit your personal needs.

The soft upholstery pampers you with thick and soft cushions – perfect for petite individuals or anyone who loves extra-comfort. Its innovative V-shape with winged structure provides optimal lumbar support, promotes relaxation, and ensures a comfortable posture.







### Effortless sitting down and getting up – comfort and safety combined

This armchair actively supports you when sitting down and standing up – gently and effortlessly at the push of a button. When raised, you can lean against it carefully as the armchair lowers you smoothly into a seated position. To stand up, simply use the remote control to gradually lift yourself into an upright position, allowing you to exit the chair with ease. The integrated lift-assist function helps prevent falls and supports your independence.

### Maximum comfort: infinitely adjustable from “feet up” to nap time

The cosy Saljol Riser Armchair is your personal oasis of relaxation at home. Whether you want to put your feet up or enjoy a restful nap, simply use the remote control to extend the footrest and recline the backrest smoothly. The independent adjustment of the footrest and backrest allows you to customise your seating and lying position to your preference. For even greater comfort, the adjustable headrest can be tilted forward or backward – just the way you like it.



### Clever & practical: the magnetic remote control

Your remote is always right where you need it! Thanks to its magnetic function, it securely attaches to the headrest or armrest of your chair – no more searching or misplacing it. With seven pre-programmed buttons, you can adjust the backrest and footrest either individually or in sync, exactly to your preference. Simple and intuitive!







### Ultimate comfort: additional accessories for your armchair

Enhance your relaxation and personalise your comfort with these three accessory highlights:

Head & lumbar cushion – the head cushion promotes a relaxed posture, while the lumbar cushion relieves the lower back and supports an ergonomic sitting position.

Tablet/smartphone holder – keep your devices in view and within easy reach – without the hassle of holding them. Enjoy movies, music, or effortless communication with family and friends from the comfort of your armchair.

Swivel table „Tabula“ – a practical and flexible table for drinks, plates, or books. Its convenient swivel function adapts to your needs.



### Always charged – USB charging station directly on the armchair

No more searching for power sockets! With two integrated USB Type-A ports, you can conveniently charge your phone or tablet right from your armchair. Whether you're making calls, chatting, or reading, enjoy your relaxation time without worrying about the battery level.



### Everything within reach in the practical side pocket

Thoughtfully designed down to the smallest detail, the integrated side pocket offers space for your most essential items. Whether it's a book, reading glasses, or the remote control – everything is safely stored and always within arm's reach. This way, you can enjoy relaxing hours without having to search for anything.







German engineering from Bavaria.

German Design



### Great luxury that remains affordable


Our riser armchairs combine the highest quality with German design – delivering maximum comfort, relaxing moments, effortless standing up, and an overall pleasant seating experience.



### Comprehensive care and support – for your mobility and safety

Our physiotherapist, Katharina, supports you in all areas where expert physiotherapeutic knowledge is required. Our helpful advisory service is always at your disposal – whether at the medical supply store or directly with us. Your safety and mobility are our top priority! If you use a rollator, our „Rollator Mobility“ training provides valuable tips and tricks for daily use. Visit us at [www.rollatormobility.com](http://www.rollatormobility.com) for more safety and enjoyment on the go!


Scan the QR codes to explore our rollator product brochures and get inspired by how a rollator can make your daily life more comfortable and secure.



Saljol offers not only high-quality riser armchairs but is also a specialist in premium rollators. These aids enhance your quality of life in everyday activities. We are committed to expanding your freedom and improving your mobility, as strong and flexible muscles are not just the key to relaxed sitting and effortless standing. To further boost your mobility, we offer targeted exercise routines to counteract prolonged sitting. Enhance your well-being with Rollastic, simple exercises that you can even perform directly in your armchair.

### Peter: „In the reclining position, I can relax and daydream!“

„I find the riser armchair incredibly comfortable. It's perfect for unwinding, and later I can move around the apartment with my rollator. The lift-assist function is especially important to me. I can raise the chair, and I'm already at a height that allows me to walk. Another key feature is that I can recline the chair to fully relax, close my eyes, and drift off into a daydream – it's just great! Plus, the chair is very easy to use.“



Saljol Riser Armchair	Firm upholstery	Soft upholstery
Seat width	47 cm	42 cm
Seat height	50 cm	50 cm
Seat depth	55.9/61 cm	55.3/60.4 cm
Recommended body height	160–190 cm	160–190 cm
Chair length in resting position	190.5 cm	190.5 cm
Maximum load	150 kg	150 kg
Space requirement	200 x 80 cm	200 x 80 cm
Armchair height	122 cm	122 cm
Armchair width	77.5 cm	77.5 cm
Inclination backrest	160°	162°

### Riser Armchair accessories



Tablet holder  
(without tablet)



Head pillow\*



Lumbar  
pillow\*



Swivel table



Incontinence  
pad set\*

\*Available in Platinum Grey and Graphite Grey

### Colour variants



Platinum Grey



Graphite Grey

Products may differ from illustrations | 07/2025