



 **SALJOL**
Joy of life!

RA
Riser Armchair

For more joy of life – premium quality with many extras

Includes remote control, side pocket and USB port

Smart and elegant
design – two
upholstery options
and colours

Quality that
impresses –
high-quality
materials for
ideal support and
pleasant softness

Maximum comfort
– continuously
adjustable from ‘feet
up’ to midday nap



Standing up and
sitting down made
easy – comfort and
safety in one

Smart and practical
– magnetic remote
control

Always charged –
USB charging station
directly on the chair

Everything within easy
reach – with side pocket

Saljol Riser Armchair	Firm upholstery	Soft upholstery
Seat width	47 cm	42 cm
Seat height	50 cm	50 cm
Seat depth	55.9/61 cm	55.3/60.4 cm
Recommended body height	160–190 cm	160–190 cm
Chair length in resting position	190.5 cm	190.5 cm
Maximum load	150 kg	150 kg
Space requirement	200 x 80 cm	200 x 80 cm
Armchair height	122 cm	122 cm
Armchair width	77.5 cm	77.5 cm
Inclination backrest	160°	162°

Riser Armchair accessories



Tablet holder
(without tablet)



Head pillow*



Lumbar
pillow*



Swivel table



Incontinence
pad set*

*Available in Platinum Grey and Graphite Grey

Colour variants



Platinum Grey



Graphite Grey

Products may differ from illustrations | 07/2025