



Saljol Riser Armchair

Comfort and safety
for your home



Highlights



One armchair – two upholstery options

Firm: comfortable seating
Soft: gentle sinking into the chair
Choose your seating comfort.
All upholstery parts can be freely
combined, are removable and
washable at 60°C. You can also
adjust the seat width without
changing the width of the chair.

Magnetic remote control
adheres to the side of the arm-
rest and the side of the headrest
– no more bending down
to reach the remote control

Daybed function
a nap in the chair is possible,
the side panel prevents you from
falling out

Adjustable seat depth
individually adjustable to suit
your leg length


Two USB ports
easily charge your phone
via the chair's built-in ports



Technical data

Saljol Riser Armchair	Firm upholstery	Soft upholstery
Seat width	47 cm	42 cm
Seat height	50 cm	50 cm
Seat depth	55.9/61 cm	55.3/60.4 cm
Recommended body height	160–190 cm	160–190 cm
Chair length in resting position	190.5 cm	190.5 cm
Maximum load	150 kg	150 kg
Space requirement	200 x 80 cm	200 x 80 cm
Armchair height	122 cm	122 cm
Armchair width	77.5 cm	77.5 cm
Inclination backrest	160°	162°

Colour variants



Platinum grey



Graphite grey

Accessories



Lumbar pillow*



Head pillow*



Incontinence pad set*



Tablet holder



Swivel table

*Available in Platinum Grey and Graphite Grey

Storage and support in one – the swivel table



Products may differ from illustrations

10/2024