

WR

Indoor Rollator



For more joy of life – premium quality and fully equipped

Includes wooden seat board and tray with felt cover, cork handrail padding, small bag and brake bar in frame colour

Stable like a mobile handrail – optimum support with every movement

Safety you can feel – skin-friendly handrail pads for maximum comfort

Domestic help for every situation – transport made easy



Twice the safety – with two braking options

Rolling seat – allows effortless movement while seated by shuffling your feet

Highly agile even in confined spaces – walk, stand and turn on the spot

Indoor Rollator	WR43	WR50	WR53	WR56
Recommended body height	135–155 cm	145–165 cm	150–180 cm	170–200 cm
Maximum user weight	120 kg	150 kg	150 kg	150 kg
Handle height	78 cm	80 cm	86 cm	92 cm
Seat height	43 cm	50 cm	53 cm	56 cm
Width bottom (wheels)	56 cm	56 cm	56 cm	56 cm
Width top	51 cm	51 cm	51 cm	51 cm
Distance between handles	43 cm	43 cm	43 cm	43 cm
Weight without accessories	9.4 kg	9.6 kg	9.8 kg	10.0 kg

Full equipment components



Small bag



Tray



Back roll



Soft seat



JOSt bag
transparent

Frame colour



Champagne Gold



Velvet Grey

Fully equipped WR: includes wooden seat and wooden tray with felt pad, handrail pad set, small bag and brake bar in frame colour

Products may differ from illustrations | 09/2025