



WR

Indoor Rollator



For more joy of life – premium quality and fully equipped

Includes wooden seat board and tray with felt cover, cork handrail padding, small bag and brake bar in frame colour

Stable like a mobile handrail – optimum support with every movement

Safety you can feel – skin-friendly handrail pads for maximum comfort

Domestic help for every situation – transport made easy



Twice the safety – with two braking options

Rolling seat – allows effortless movement while seated by shuffling your feet

Highly agile even in confined spaces – walk, stand and turn on the spot

Indoor Rollator	WR43	WR50	WR53	WR56	WRXL*
Recommended body height in cm	135–155	145–165	150–180	170–200	170–210
Maximum user weight	120 kg	150 kg	150 kg	150 kg	180 kg
Handle height	78 cm	80 cm	86 cm	92 cm	95 cm
Seat height	43 cm	50 cm	53 cm	56 cm	58 cm
Width bottom (wheels)	56 cm	56 cm	56 cm	56 cm	68 cm
Width top	51 cm	51 cm	51 cm	51 cm	65 cm
Distance between handles	43 cm	43 cm	43 cm	43 cm	60 cm
Weight without accessories	9.4 kg	9.6 kg	9.8 kg	10.0 kg	13.4 kg
Compatible with accessories	✓	✓	✓	✓	✗

Full equipment components



Small bag



Tray



Back roll

Additional accessories



Soft seat



JOST bag
beige



JOST bag
transparent

Frame colour



Velvet Grey



Champagne Gold

*WRXL is supplied with seat board and brake bar in frame colour and is only available in Velvet Grey.