



Indoor Rollator

**Completely safe at home.
Supports without drawing attention.**

To help prevent falls
at home

At first glance: a stylish
serving trolley

As stable as a mobile
handrail, it provides
support at every point

Fits through any door
and can be turned on
the spot



Unique: move while seated
thanks to the flip-up seat

Dual brake system with
parking function:
hand brake and brake bar

Includes accessories:


- Wooden flip-up seat with felt pad
- Tray with non-slip felt pad
- Small bag
- Cork handrail pads

Indoor Rollator (WR)	WR43	WR50	WR53	WR56
Recommended body height	135–155 cm	145–165 cm	150–180 cm	170–200 cm
Maximum user weight	120 kg	150 kg	150 kg	150 kg
Handle height	78 cm	80 cm	86 cm	92 cm
Seat height	43 cm	50 cm	53 cm	56 cm
Width bottom (wheels)	56 cm	56 cm	56 cm	56 cm
Width top	51 cm	51 cm	51 cm	51 cm
Distance between handles	43 cm	43 cm	43 cm	43 cm
Weight without accessories	9.4 kg	9.6 kg	9.8 kg	10.0 kg

Accessories



Small bag



Tray with felt pad



Back roll




Soft seat

Frame colours



Champagne
Gold



Velvet Grey

Fully equipped WR: includes wooden seat and wooden tray with felt pad, handrail pad set, small bag and brake bar in frame colour