



Carbon Rollator Superlight

Probably the world's lightest premium rollator

Consistently
lightweight in every
component – only
4.4/4.5 kg

Uncompromisingly
stable and safe

Cork-TPE handles
– pleasantly soft
and comfortable
to grip

Fidlock magnetic
shopping net for
intuitive attachment



Shopping net
supports loads of up
to 10 kg with optimal
weight distribution

Exceptionally easy
to lift and store

Smooth rolling
without drawing
attention

TPE tyres with
reliable braking
performance and a
light tread pattern

Carbon Rollator (CR) Superlight	CR54 Superlight	CR62 Superlight
Recommended body height	145–170 cm	160–200 cm
Maximum user weight	150 kg	150 kg
Total weight without accessories	4.4 kg	4.5 kg
Maximum load of net bag	10 kg	10 kg
Width open/closed	61/21 cm	61/21 cm
Seat width	46 cm	46 cm
Seat height	54 cm	62 cm
Handle height (adjustable)	76–88 cm	81–103 cm

Accessories



Back belt




Firm seat



Wooden tray
with felt pad



Bag
leatherette brown



Umbrella
caramel



Umbrella
black



Grinding brake



LED light




Travel bag



Brake cable for
the conversion to
one-hand brake



Stick holder



Bell



Frame
colour

Rose Gold