



Allround Rollator

Enjoy your daily adventures

100% stable on any terrain: only 6.1/6.3 kg lifting weight

Very easy to fold, compact to store

Double brake package:

- Hand brake (convertible to a one-hand brake)
- Parking brake

All-inclusive package:

- Bell
- Stick holder



Comfort package:

- Premium shoulder bag with magnetic flap
- Fidlock magnetic attachment for intuitive fastening and 10 kg load capacity
- Back belt with cut-out
- Ergonomic, warm cork-TPE handles
- Large, puncture-proof all-terrain tyres for smooth obstacle clearance and safe manoeuvring, even on uneven paths

Visibility package:

- 30 reflectors
- LED light

Allround Rollator (AR)	AR54	AR62	AR62 XL*
Recommended body height	145–170 cm	160–200 cm	160–200 cm
Maximum user weight	150 kg	150 kg	180 kg
Total weight without accessories	6.1 kg	6.3 kg	6.4 kg
Maximum load of bag	10 kg	10 kg	10 kg
Width open/closed	61/21 cm	61/21 cm	68/21 cm
Seat width	46 cm	46 cm	52 cm
Seat height	54 cm	62 cm	62 cm
Handle height (adjustable)	76–88 cm	81–103 cm	81–100 cm

*AR62 XL only available in Barolo Red

Accessories



Back belt



Firm seat



Wooden tray
with felt pad



Net bag
grey



Umbrella
caramel



Umbrella
black



Grinding brake



LED light




Travel bag



Brake cable for
the conversion to
one-hand brake



Stick holder



Bell

Frame colours



Pearl Grey



Magic Black



Forest Green



Barolo Red