



Carbon Rollator

Elegance and freedom on four wheels

Pleasantly lightweight: only
5.3/5.5 kg lifting weight

Very easy to fold,
compact to store

Triple brake package:

- Hand brake (convertible to a one-hand brake)
- Parking brake
- Grinding brake for safe walking downhill and to counteract a left or right drift

All-inclusive package:

- Bell
- Stick holder



Comfort package:

- Premium shoulder bag with magnetic flap
- Fidlock magnetic attachment for intuitive fastening and 10 kg load capacity
- Back belt with cut-out
- Ergonomic, warm cork-TPE handles
- Soft, puncture-proof tyres absorb shocks and provide secure grip

Visibility package:

- 30 reflectors
- LED light

Carbon Rollator (CR)	CR54	CR62	CR62 XL*
Recommended body height	145–170 cm	160–200 cm	160–200 cm
Maximum user weight	150 kg	150 kg	180 kg
Total weight without accessories	5.3 kg	5.5 kg	5.6 kg
Maximum load of bag	10 kg	10 kg	10 kg
Width open/closed	61/21 cm	61/21 cm	68/21 cm
Seat width	46 cm	46 cm	52 cm
Seat height	54 cm	62 cm	62 cm
Handle height (adjustable)	76–88 cm	81–103 cm	81–100 cm

*CR62 XL only available in British Racing Green

Accessories



Back belt



Firm seat



Wooden tray
with felt pad



Net bag
grey



Umbrella
caramel



Umbrella
black



Travel bag



LED light



Stick holder



Bell



Brake cable for
the conversion to
one-hand brake

Frame colours



British Racing
Green



Midnight Blue



Copper Red



Pure Gold