



Allround Rollator Superlight

Brings lightness to your adventures

Extremely stable yet exceptionally lightweight in its class – only 5.2/5.3 kg

Large front wheels – ideal even for uneven terrain

Cork-TPE handles – pleasantly soft and comfortable to grip

For effortless everyday adventures

Fidlock magnetic shopping net for intuitive attachment



Shopping net supports loads of up to 10 kg with optimal weight distribution

Enjoys an adventure too

Easy to fold and convenient to transport

TPE tyres with reliable braking performance and a light tread pattern

An allround rollator with character and sophistication

Allround Rollator (AR) Superlight	AR54 Superlight	AR62 Superlight
Recommended body height	145–170 cm	160–200 cm
Maximum user weight	150 kg	150 kg
Total weight without accessories	5.2 kg	5.3 kg
Maximum load of net bag	10 kg	10 kg
Width open/closed	61/21 cm	61/21 cm
Seat width	46 cm	46 cm
Seat height	54 cm	62 cm
Handle height (adjustable)	76–88 cm	81–103 cm

Accessories



Back belt




Firm seat



Wooden tray
with felt pad



Bag
leatherette brown



Umbrella
caramel



Umbrella
black



Grinding brake



LED light




Travel bag



Brake cable for
the conversion to
one-hand brake



Stick holder



Bell

Frame
colour



Ice Blue