



Shower Chair

Safe and comfortable all-round showering experience

Adjustable seat and back unit with braking function when leaning on the armrests

360° rotatable for greater freedom of movement and safe turning even in small showers

Secure grip thanks to 5 anti-slip castors for controlled rotational movement

Personal care without getting wet yourself, thanks to smooth 360° rotation



Heat-reflecting, adjustable seat and back surfaces for comfortable personal care thanks to hygiene cut-out

Stable yet lightweight construction made of rust-free aluminium: only 5.5 kg, supports up to 150 kg

Versatile use in the bathroom, also suitable as a hairdressing or make-up chair

Shower Chair (DS)	DS50
Recommended body height	150–190 cm
Maximum user weight	150 kg
Seat height	50 cm
Hygiene cut-out, position 1/2/3	6/9/12 cm
Seat width	36/39/42 cm
Seat depth	37 cm
Backrest length	37 cm
Backrest height	33/36/39 cm
Required turning space (product width)	60 cm
Total weight without accessories	5.5 kg

Accessories

Adapter set for seat height increase (+4/+8 cm) and seat angle adjustment (13°)



Adapter set for adjusting the seat height to 54 or 58 cm, or for tilting the seat surface to create an ergonomic wedge seating position



Frame colour



Diamond Silver