



Figure: Shower Stool with backrest

Shower Stool

**Showering experience like standing.
Turn on the spot while seated.**

360° rotatable for greater freedom of movement and safe turning even in small showers

Secure grip thanks to 5 anti-slip castors for controlled rotational movement

Personal care without getting wet yourself, thanks to smooth 360° rotation



Heat-reflecting, adjustable seat surfaces for comfortable personal care thanks to hygiene cut-out

Stable yet lightweight construction made of rust-free aluminium: only 2.8 kg, supports up to 150 kg

Versatile use in the bathroom, also suitable as a hairdressing or make-up stool

Shower Stool (DH)	DH50
Recommended body height	150–190 cm
Maximum user weight	150 kg
Seat height	50 cm
Hygiene cut-out, position 1/2/3	6/9/12 cm
Seat width	36/39/42 cm
Seat depth	37 cm
Required turning space (product width)	52 cm
Total weight without accessories	2.8 kg

Accessories

Backrest



Backrest to support a safe and upright seating position



Adapter set for seat height increase (+4/+8 cm) and seat angle adjustment (13°)



Adapter set for adjusting the seat height to 54 or 58 cm, or for tilting the seat surface to create an ergonomic wedge seating position

Frame colour



Diamond Silver